



## HISTOPRIME<sup>®</sup>

**CatNo E131**

### TBS-Buffer pH 7,2

**Lot: See Label**  
**Storage: +2 to +8°C**  
**Exp. Date: See Label**

**Composition**

50 mM buffer salts of tris-(hydroxymethyl)-aminomethane (base) and HCL, 015 M NaCl

**pH value**

7.2 at 20°C

**Content**

500 ml ready-to-use buffers

**Usage**

Can be used as wash buffer for staining with HistoGreen (LINARIS Cat. No. E109).

**Storage**

In the refrigerator at 2 -8°C.

E131 230109-1/1

**For Research purposes only. Not for therapeutic or diagnostic use.**

