



## Normal Donor Horse Serum, sterile

**CatNo.:** ENHD2000

**Lot Number:** See Label

**Size:** CatNo.: ENHD2000-50; 50ml  
CatNo.: ENHD2000-500; 500ml

**Application:** Cell Culture  
Immunohistology – Frozen  
Immunohistology – Paraffin  
Immunofluorescence  
Western Blotting  
Blocking  
Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures.

**Purity:** Whole antiserum. Final filter size 0.2µm.

**General notes:** Normal Donor Horse serum can be used for blocking non-specific conserved sequence and/or Fc receptor binding.  
Typically, the serum used for blocking is from a different species than the species in which the primary antibody was raised. Often the blocking serum is from the species in which the secondary antibody was raised.

**Collection:** When searchers choose their serum, an important factor that should be taken into consideration is the source, which also emphasizes the traceability of the serum. Our system of vertical integration allows us to be certain of the origin and traceability of our serum. The serum is imported and treated in agreement with the European regulations. Mixed breed and sex.

**Virus Tested:** All of our serum is tested for the Equine Infectious Anemia (EIA)

**Sterility:** All sera are tested for the absence of aerobic and anaerobic bacteria, fungi, yeast and Mycoplasma. The sterility test is based on the European Pharmacopeia requirements.  
The sera are tested for the absence of Mycoplasma by culture.

**Protein Concentration:** 74.4 mg/ml

**Shipping:** Ship on Dry Ice.

**Storage:** -20°C

**Shelf Life:** 12 months from date of dispatch.

**Regulatory:** For **in vitro research** use only, not for diagnostic or therapeutic use.

Manufacturer

**BIOZOL**  
FIT FOR SCIENCE

**BIOZOL**

Diagnostica Vertrieb GmbH  
Leipziger Straße 4  
85386 Eching

Phone +49 (89) 3799 666 6  
Fax +49 (89) 3799 666 99  
E-Mail info@biozol.de  
www biozol.de



Page 1 of 1