

## HISTOPRIME<sup>®</sup>

CatNo LHS0202

### HistoClear II (XEM-200)

Lot: See Label  
Storage: 24°C  
Exp. Date: See Label

**Solvent substitute (xylene substitute) and clarifying agent based on citrus oil.**

<b>General</b>	Non-toxic solution for laboratory use, with citrus aroma. XEM-200 can be used as a substitute for xylene and toluene. It can be used in dyeing lines and in embedding machines. XEM-200 is suitable for dissolving kerosene and fats. Therefore, to prevent drying of the skin, direct contact with XEM-200 should be avoided.
<b>Packing</b>	3785 ml (1 Gal)
<b>Components</b>	100% natural citrus oil extracts. Without chemical additives.
<b>Storage</b>	Store carefully closed, in a cool and dry environment.
<b>Precautions</b>	<b>Not suitable for human consumption.</b>

LHS0202 210901 1/1



Reizend

Hinweise: R 10-38  
S 28  
XEM-200 an Händen oder in den Augen ist  
sofort mit klarem Wasser zu entfernen.  
Flammpunkt: 49°C

**For Research purposes only. Not for therapeutic or diagnostic use.**