

## Bovine Serum Albumin (BSA),

I

(IgG-Free, Protease-Free)

CatNo

E140

BatchNo: See Label Storage: 2 - 25°C

Amount:	50 g
Remarks:	This highly purified grade of bovine serum albumin is used in all of our products which contain BSA as a carrier. Our IgG-Free and Protease -Free grade of bovine serum albumin is recommended for use in diluents for antibodies and enzymes, or as a blocking agent for many immunhistochemical and transfer blot applications.
Purity:	Tested negative for bovine IgG and for protease.
Expiration Date:	Three Years from date of receipt. The expiration date may be extended if test results are acceptable for the intended use.
Rehydration:	Weigh out desired amount and dissolve in desired buffer. Centrifuge if not clear. Prepare working dilution on day of use.
Storage:	2 - 25°C in dry conditions protected from light.
E140 231211-1/1	This product is for in vitro research use only. It is not a medical device and it I not in- tended for diagnostic or therapeutic purposes.



BIOZOL Diagnostica Vertrieb GmbH Leipziger Straße 4 85386 Eching

Phone +49 (89) 3799 666 6 +49 (89) 3799 666 99 E-Mail info@biozol.de biozol.de

Fax

www



