



Bovine Serum Albumin 30%, sterile

CatNo LIN0297-100

Lot: See Label
Storage: 2-8°C
Exp. Date: See Label

Amount: 100 ml
Solution: 30% BSA, the solution is prepared and filtered through a 0.2 µm filter.
Shelf life: 12 months
Sterility tests:
- Bacteria in aerobic and anaerobic conditions
- Fungi and yeasts

Recommended use: Wide range of applications, e.g. as stabilizing proteins in solution, blocking non-specific binding in immuno-assays, and as a protein standard in various protein assays.
BSA 30% solution is a convenient form for direct supplementation of serum-free media. Bovine serum albumin is broadly used as an additive to cell culture media, especially serum-free media. It provides a range of benefits including protection from oxidative damage and stabilization of other media components such as fatty acids and pyridoxal.

Storage condition: Store at 4°C. Stable for at least 12 months from date of receipt when stored as recommended.

LIN0297-100-230110-1/1

For Research Use Only

Manufacturer

BIOZOL
FIT FOR SCIENCE

BIOZOL

Diagnostica Vertrieb GmbH
Leipziger Straße 4
85386 Eching

Phone +49 (89) 3799 666 6
Fax +49 (89) 3799 666 99
E-Mail info@biozol.de
www biozol.de



Page 1 of 1